



Malosa Medical
Single-Use Surgical Instruments

Aesthetic Surgery

www.malosa.com

ISO Certificates

**All Malosa instruments and procedure packs are
cleaned, packed and sterilised in the UK
under UK accredited ISO standards**



ISO 13485:2003



ISO 9001:2000



**93/42 EEC
CE Marking**



Malosa Medical

SINGLE-USE INSTRUMENTS

CONTENTS

ii-iv **Single Use Aesthetic Instruments List**

INSTRUMENTATION

Breast Retractor

Abdominal Liposuction Cannula

Facial Liposuction Cannula

Infiltration and Fat Repositioning Cannula

Scissors

Ramplays

Forceps

Retractors and Skin Hooks

Haemostatic Forceps

Needle Holders

MALOSA PACKS

Breast Augmentation Pack

Rhytidectomy Pack

Blepharoplasty Pack

Liposuction Pack

Ocular Plastics Pack

Chalazion Pack

Minor Plastics Pack

Minor Ops Hand Pack

Suturing and Suture Removal Packs



Single-Use Aesthetic Instruments List

BREAST RETRACTORS

- MMSU3012 Breast Retractor with Fibre Optics and Smoke Aspiration Tube, 25mm blade
- MMSU3013 Breast Retractor with Fibre Optics and Smoke Aspiration Tube, 25mm toothed blade
- MMSU3014 Breast Retractor with Fibre Optics and Smoke Aspiration Tube, 15mm blade

LIPOSUCTION CANNULA

Infiltration

- MMSU3217 Liposuction Infiltration Cannula, 2.75mm ID, 325mm Long, Standard tip, 8 Spiral ports, Luer Lock
- MMSU3218 Liposuction Infiltration Cannula, 3.0mm ID, 200mm Long, Standard tip, 4 Spiral ports, Luer Lock
- MMSU3219 Liposuction Infiltration Cannula, 2.0mm ID, 185mm Long, Standard tip, 8 Spiral ports, Luer Lock
- MMSU3230 Liposuction Infiltration Cannula, 2.0mm OD 100mm Long, 8 Spiral ports, Luer Lock
- MMSU3231 Liposuction Infiltration Cannula 2.0mm OD, 150mm Long, 8 Spiral ports, Luer Lock
- MMSU3232 Liposuction Infiltration Cannula 2.0mm OD, 200mm Long, 8 Spiral ports, Luer Lock
- MMSU3233 Liposuction Infiltration Cannula 2.0mm OD, 250mm Long, 8 Spiral ports, Luer Lock

Abdominal

- MMSU3103 Liposuction Cannula, 3.0mm Dia x 230mm Long, Standard tip, Tri-port Radial
- MMSU3104 Liposuction Cannula, 3mm Dia X 230mm, Standard Tip, Pyramid Port
- MMSU3105 Liposuction Cannula, 4.0mm Dia x 260mm Long, Standard tip, Tri-Port Radial
- MMSU3106 Liposuction Cannula, 4.0mm Dia x 260mm Long, Standard tip, Pyramid Port
- MMSU3107 Liposuction Cannula, 5.0mm Dia x 260mm Long, Standard tip, Pyramid Port - Luer Lock
- MMSU3108 Liposuction Cannula, 5mm Dia X 260mm, Standard Tip, Tri-Port Radial
- MMSU3113 Liposuction Cannula, 6mm Dia X 275mm, Standard Tip, Pyramid Port
- MMSU3125 Liposuction Cannula, 4.0mm Dia x 260mm Long, Bullet tip, Tri-port radial
- MMSU3126 Liposuction Cannula - 3.0mm Dia, 230mm, Bullet Tip, Tri-Port Radial, 5.5 x 2.0mm ports
- MMSU3127 Liposuction Cannula, 5mm Dia X 260mm, Bullet Tip, Tri-Port Radial
- MMSU3143 Liposuction Cannula, 3mm Dia X 260mm, Spatula Tip, Mono Port
- MMSU3145 Liposuction Cannula, 4mm Dia X 260mm, Spatula Tip, Mono Port

Facial

- MMSU3100 Titanium Luer Lock F/F Connectors
- MMSU3119 Liposuction Cannula, 2.0mm ID, 70mm Long, Standard tip, Tri-Port Radial, Luer lock
- MMSU3120 Liposuction Cannula, 3.0mm ID, 110mm Long, Standard tip, Tri-Port Radial, Luer lock
- MMSU3122 Liposuction Cannula, 2.5mm ID, 110mm Long, Standard tip, Tri-Port Radial, Luer lock
- MMSU3123 Liposuction Cannula, 2.75mm ID, 185mm Long, Standard tip, Tri-Port Radial, Luer lock
- MMSU3124 Liposuction Cannula, 2.75mm ID, 185mm Long, Standard tip, Bi-Port Radial, Luer lock
- MMSU3141 Liposuction Cannula, 2.5mm ID, 185mm Long, Spatula tip, Mono Port, Luer Lock
- MMSU3142 Liposuction Cannula, 2.5mm ID, 135mm Long, Spatula tip, Mono Port, Luer Lock
- MMSU3159 Liposuction Cannula, 2mm X 80mm, Standard Tip, Mono Port, Luer Lock
- MMSU3165 Liposuction Cannula, 2.0mmID, 100mm long, spatula tip, Bi Port Axial, Luer Lock
- MMSU3166 Liposuction Cannula, 2.0mmID, 150mm long, spatula tip, Bi Port Axial, Luer Lock
- MMSU3167 Liposuction Cannula, 2.0mmID, 200mm long, spatula tip, Bi Port Axial, Luer Lock
- MMSU3168 Liposuction Cannula, 2.0mmID, 250mm long, spatula tip, Bi Port Axial, Luer Lock
- MMSU3171 Liposuction Cannula, 2.0mm ID, 135mm Long, Spatula tip, Bi-Port Axial, Luer Lock
- MMSU3172 Liposuction Cannula, 2.0mm ID, 185mm Long, Spatula tip, Bi-Port Axial, Luer Lock
- MMSU3173 Liposuction Cannula, 3.0mm ID, 135mm Long, Spatula tip, Bi-Port Axial, Luer Lock
- MMSU3174 Liposuction Cannula, 3.0mm ID, 185mm Long, Spatula tip, Bi-Port Axial, Luer Lock
- MMSU3179 Liposuction Cannula, 2mm X 150mm, Standard Tip, Bi-Port Axial, Luer lock
- MMSU3181 Liposuction Cannula, 3.0mm Dia x 150mm Long, Standard tip, Bi-port Axial, Luer Lock
- MMSU3199 Liposuction Cannula, 3mm X 80mm, Standard Tip, Tri-Port Axial, Luer Lock
- MMSU3220 7 Piece Facial Liposuction Cannula Set, Titanium Luer Lock 210 mm X 3.0mm - Bi Port Radial
- MMSU3225 Liposuction Fat Repositioning Cannula, 2.0mm Dia, 80mm Blunt Tip, 4mm Port, Luer Lock
- MMSU3226 Liposuction Fat Repositioning Cannula, 2.0mm Dia, 90mm Blunt Tip, 4mm Port, Luer Lock



Single-Use Aesthetic Instruments List

SCISSORS

Metzenbaum

- MMSU0243 Curved Metzenbaum Scissors, 11.5cm
- MMSU0242 Curved Metzenbaum Scissors, 18cm
- MMSU0254 Straight Metzenbaum Scissors, 18cm
- MMSU0260 Metzenbaum Rhytidectomy Scissors, 18cm
- MMSU0261 Metzenbaum Rhytidectomy Scissors, 23cm
- MMSU0262 Metzenbaum Rhytidectomy Scissors, 26cm

Stephens Tenotomy

- MMSU0208 Stephens Tenotomy Scissors, Curved, Blunt tip
- MMSU0209 Stephens Tenotomy Scissors, Curved, Sharp tip
- MMSU0213 Stephens Tenotomy Scissors, Straight, Sharp tip
- MMSU0214 Stephens Tenotomy Scissors, Straight, Blunt tip

Iris Dissecting

- MMSU0244 Iris Dissecting Scissors, Curved, Sharp tip
- MMSU0249 Iris Dissecting Scissors, Straight, Sharp tip
- MMSU0294 Iris Dissecting Scissors, Curved, Blunt tip
- MMSU0299 Iris Dissecting Scissors, Straight, Blunt tip
- MMSU0268 Iris Dissecting Scissors, Straight, Sharp tip - Gold handle
- MMSU0270 Iris Dissecting Scissors, Curved, Sharp tip - Gold handle

General

- MMSU0210 Cottle Nasal Cartilage Scissors, 5cm blade- Gold handle
- MMSU0215 Graefe Scissors - Gold handle
- MMSU0211 Kilner Scissors, Curved, Blunt tip, 16.5cm
- MMSU0256 Reynolds Dissecting Scissors, Curved, 15cm
- MMSU0255 Mayo Scissors 14.5cm
- MMSU0218 Standard Scissors
- MMSU0153 Large Blunt Scissors
- MMSU0140 Curved Blade Scissors, Blunt tip, 9cm
- MMSU0141 Curved Blade Scissors, Sharp tip, 9cm
- MMSU0154 Straight Dissecting Scissors, Sharp tip, 9cm
- MMSU0155 Straight Dissecting Scissors, Blunt tip, 9cm

FORCEPS

Non-Toothed

- MMSU0080 Fine Toothed Forceps, 10.4cm
- MMSU0278 Iris Toothed Forceps, 10.5cm
- MMSU0238 Adson Tissue Forceps, 1.25mm tip, 12.5cm
- MMSU0231 Waughs Tissue Forceps, 18cm
- MMSU0279 St Martins Forceps, 8.5cm
- MMSU0276 Bishops Forceps, 8.5cm
- MMSU3010 Converse Green Fixation Forceps

Non-Toothed

- MMSU0247 Iris Forceps, 11.5cm
- MMSU0241 Adson Dissecting Forceps, 1.25mm tip, 12.5cm
- MMSU0230 Waughs Dissecting Forceps, 18cm
- MMSU0267 Iris Dissecting Forceps, Half Curved, Serrated
- MMSU0277 Moorfields Forceps, 10.5cm
- MMSU0085 Flat Tipped Forceps, 8cm



Single-Use Aesthetic Instruments List

RAMPLEYS

- MMSU0068 Plastic Holding Rampley, Straight, 18cm
- MMSU0069 Plastic Holding Rampley, Angled tip, 11cm
- MMSU0050 Rampley Sponge Holding Forceps, 25cm
- MMSU3014 Backhaus Towel Clip
- MMSU0089 Jones Towel Clamp, 9cm
- MMSU0199 Bulldog Clip

RETRACTORS AND SKIN HOOKS

- MMSU3001 Desmarres Lid Retractor
- MMSU3031 Lagenbach Retractor, 21.5cm handle
- MMSU3035 Czerny Retractor, 16.5cm
- MMSU0252 Muscle Hook
- MMSU3003 West Retractor 3 x 4 teeth, Sharp
- MMSU3005 Guthrie Hook, Fine Hooks, 2mm wide
- MMSU3006 Gillies Skin Hook, 13cm
- MMSU3037 Gillies Skin Hook, 16.5cm
- MMSU3036 McIndoe Skin Retractor, 10mm Sharp, 15cm Long

HAEMOSTATIC FORCEPS

- MMSU3007 Halstead Mosquito Forceps, Curved
- MMSU3008 Halstead Mosquito Forceps, Straight
- MMSU0077 Spencer Wells Forceps, Straight, 23cm
- MMSU0051 Dunhill Forceps, 13cm

FORCEPS

- MMSU0182 Crile-Murray Needle Holder, without groove, 13cm
- MMSU0183 Derf Needle Holder, without groove, 12cm
- MMSU0184 Crile-Murray Needle Holder, without groove, 15cm
- MMSU0185 Halsey Needle Holder
- MMSU3012 Castroviejo Needle Holder, Straight with lock
- MMSU3013 Castroviejo Needle Holder, Curved with lock



Malosa Medical
Single-Use Surgical Instruments

Instrumentation

www.malosa.com




Malosa Medical

SINGLE-USE INSTRUMENTS

Fibre Optic Retractor

All surgical steel retractor with genuine fibre optics and aspiration tube

3 Options



3412 - 25mm Wide



3413 - 25mm Wide with Teeth



3414 - 15mm Wide



Rigid lightweight construction

Re-usable quality at disposable prices

New optics every time

Aspiration tube



Abdominal Liposuction Cannula

Available in a range of tips, port configurations, lengths and diameters.



Pressure release vent on thumb grip

Bespoke designs available for a minimum order of 250pcs


Pyramid



Ref #	Length (mm)	Dia (mm)	Port (mm)
3104	230	3.0	5.5 X 2.3
3106	260	4.0	7.0 X 2.5
3113	275	6.0	8.5 X 3.0

Mercedes

Standard Tip



Ref #	Length (mm)	Dia (mm)	Port (mm)
3103	230	3.0	5.5 X 2.0
3105	260	4.0	7.0 X 2.25
3108	260	5.0	8.5 X 2.5

Bullet Tip



Ref #	Length (mm)	Dia (mm)	Port (mm)
3126	230	3.0	5.5 X 2.0
3125	260	4.0	7.0 X 2.0
3127	260	5.0	8.5 X 2.3



Facial Liposuction Cannula

Cannula with leur lock hub, ports in a radial orientation, with a range of tips, port numbers, lengths and diameters.



Standard Tip

Mono-Port



Ref #	Length (mm)	Dia (mm)	Port (mm)
3159	80	2.0	3.5 X 1.0

Bi-Port Radial



Ref #	Length (mm)	Dia (mm)	Port (mm)
3124	185	2.75	4.0 X 2.0

Tri-Port Radial



Ref #	Length (mm)	Dia (mm)	Port (mm)
3119	70	2.0	3.0 X 1.0
3120	110	3.0	4.0 X 1.5
3122	110	2.5	3.0 X 0.8
3123	185	2.75	4.0 X 2.0



Facial Liposuction Cannula

Cannula with leur lock hub, ports aligned on the axis with a range of tips, port numbers, lengths and diameters.


Standard Tip

Bi-Port Axial



Ref #	Length (mm)	Dia (mm)	Port (mm)
3179	150	2.0	3.0 X 1.5
3181	150	3.0	6.0 X 2.0

Tri-Port Axial



Ref #	Length (mm)	Dia (mm)	Port (mm)
3199	80	3.0	3.5 X 2.0

Spatula Tip

Mono-Port



Ref #	Length (mm)	Dia (mm)	Port (mm)
3142	135	2.5	3.0 X 1.5

Bi-Port Axial



Ref #	Length (mm)	Dia (mm)	Port (mm)
3171	135	2.0	5.0 X 1.0
3172	185	2.0	5.0 X 1.0
3173	135	3.0	6.5 X 2.0
3174	185	3.0	6.5 X 2.0



Liposuction Cannula

Cannula with leur lock for tumescent infiltration and fat repositioning

Fat Repositioning Cannula

Extra wide and longer port to prevent fat cell damage



Ref #	Length (mm)	Dia (mm)	Port (mm)
3225	80	2.0	4.0 X 1.5
3226	90	2.0	4.0 X 1.5

Infiltration Cannula

4 Spiral Port




Ref #	Length (mm)	Dia (mm)	Port (mm)
3218	200	3.0	1.0 Dia

8 Spiral port



Ref #	Length (mm)	Dia (mm)	Port (mm)
3217	325	2.75	1.0 Dia
3219	185	2.0	1.0 Dia

Leur Lock Connector



Ref #	Length (mm)	Dia (mm)	Port (mm)
3100	20	9.0	4.0 Dia



Malosa Medical

SINGLE-USE INSTRUMENTS

Scissors

Metzenbaum



Stephens Tenotomy

Iris



General





Scissors

General



Ramplays

Plastic and Steel



Forceps

Toothed Forceps



Non-Toothed Forceps





Retractors and Skin Hooks



Haemostatic Forceps



Needle Holders





Malosa Medical
Single-Use Surgical Instruments

Instrument Packs

www.malosa.com



Malosa Medical

SINGLE-USE INSTRUMENTS

Fibre Optic Retractor

All surgical steel retractor with genuine fibre optics and aspiration tube

3 Options



3412 - 25mm Wide



3413 - 25mm Wide with Teeth



3414 - 15mm Wide



Rigid lightweight construction

Re-usable quality at disposable prices

New optics every time

Aspiration tube



Cosmetic Surgery Packs



Breast Augmentation



MMK048

Contents:

- MMSU0260 - Metzenbaum Rhytidectomy Scissors, Curved, 18cm
- MMSU0159 - Mayo Scissors Straight, 16.5cm
- MMSU3036 - McIndoe Double Skin Hook, 16.5cm x 2
- MMSU0264 - Nelson Scissors, Curved, 22.5cm
- MMSU0218 - Standard Scissors Blunt/Blunt, 13cm
- MMSU0051 - Dunhill Artery Forceps, Curved, 13cm x 2
- MMSU3008 - Halstead Mosquito Forceps, Straight, 12.5cm x 2
- MMSU0050 - Rampley Sponge Holding Forceps, 25cm
- MMSU0238 - Adson Tissue Forceps, Lightweight, 12.5cm x 2
- MMSU0224 - Debakey Dissecting Forceps, 18cm
- MMSU0228 - Allis Tissue Forceps 3/4 Teeth, 15cm x 2
- MMSU3020 - BP No. 3 Scalpel Blade Handle with Gardations in cm
- MMSU0182 - Crile-Murray Needle Holder, 13cm
- MMSU0183 - Derf Needle Holder, 12cm x 2
- MMSU0184 - Crile-Murray Needle Holder, 15cm
- MMSU3005 - Guthrie Double Skin Hook, 16.5cm x 2
- MMSU3037 - Gillies Skin Hook, 13cm x 2
- MMSU3031 - Lagenbach Retractor, 21.5cm
- MMSU3035 - Czerny Retractor, 18cm
- MMSU0118 - Kidney Dish
- MMSU0515 - Trolley Drape, Reinforced Plastic Cover, 127 x 228cm

Instruments

 0260 Metzenbaum	 0159 Mayo	 3036 McIndoe	 0264 Nelson	 0218 Standard
 0051 Dunhill	 3008 Halstead	 0050 Rampley	 0238 Adson	 0224 Debakey
 0228 Allis	 3020 BP No. 3	 0182 Crile-Murray 13cm	 0183 Derf	 0184 Crile-Murray 18cm
 3005 Guthrie	 3037 Gillies	 3031 Lagenbach	 3035 Czerny	



Cosmetic Surgery Packs

Rhytidectomy



Contents:

- MMSU0051 - Dunhill Artery Forceps, Curved, 12.8cm
- MMSU3008 - Halstead Mosquito Forceps (Micro), Straight x 2
- MMSU0077 - Spencer Wells Artery Forceps, Straight, 23cm x 2
- MMSU3020 - No. 3 Scalpel Handle with Gradations x 2
- MMSU0230 - Waughs Dissecting Forceps, 18cm
- MMSU0231 - Waughs Tissue Forceps, 18cm
- MMSU0238 - Adson Tissue Forceps, 12.5cm
- MMSU3035 - Czerny Retractor, 18cm
- MMSU3031 - Lagenbach Retractor, 21.5cm Handle
- MMSU3036 - McIndoe Double Skin Hook, 16.5cm x 2
- MMSU3037 - Gillies Skin Hook, 16.5cm x 2
- MMSU0159 - Mayo Scissors, Straight, 16.5cm
- MMSU0211 - Kilner Scissors, Curved, 16.5cm
- MMSU0242 - Metzenbaum Scissors, Curved, 18cm
- MMSU0254 - Metzenbaum Scissors, Straight, 18cm
- MMSU0256 - Reynolds Dissecting Scissors, Curved, 15cm
- MMSU0244 - Iris Dissecting Scissors, Curved, 11.5cm
- MMSU0243 - Metzenbaum Scissors, Curved, 11.5cm
- MMSU0184 - Crile-Murray Needle Holder, 15cm
- MMSU0183 - Derf Needle Holder, 12cm
- MMSU0182 - Crile-Murray Needle Holder, 13cm
- MMSU0260 - Metzenbaum Rhytidectomy Scissors, 18cm
- MMSU0261 - Metzenbaum Rhytidectomy Scissors, 23cm
- MMSU0262 - Metzenbaum Rhytidectomy Scissors, 26cm
- MMSU0218 - Blunt Nurse's Scissors, 13cm

MMK047

Instruments



0051 Dunhill



3008 Halstead



0077 Spencer Wells



3020 BP No. 3



0230 Waughs



0231 Waughs



0238 Adson



3035 Czerny



Malosa Medical

SINGLE-USE INSTRUMENTS

Cosmetic Surgery Packs



Rhytidectomy

Instruments





Cosmetic Surgery Packs



Blepharoplasty



Contents:

- MMSU3001 - Desmarres Lid Retractor, 11mm wide
- MMSU3006 - Gillies Skin Hook, 13cm
- MMSU3005 - Guthrie Hook, 2.0mm Wide x 2
- MMSU3010 - Converse Green Fixation Forceps, 9mm Jaw
- MMSU3014 - Backhaus Towel Clip, 11cm
- MMSU3012 - Castroviejo Needle Holder - with Lock, 9mm Jaw
- MMSU3020 - BP No.3 Scalpel Handle with Gradations in CM
- MMSU1124 - Fine Toothed Forceps, 0.4mm 1 x 2 Teeth x 2
- MMSU1401 - Castroviejo Caliper, to 20mm
- MMSU3007 - Halstead Mosquito Forceps, Curved x 2
- MMSU0244 - Iris Dissecting Scissors, Curved, 11.5cm
- MMSU0249 - Iris Dissecting Scissors, Straight, 11.5cm
- MMSU3008 - Halstead Mosquito Forceps, Straight
- MMSU0215 - Graefe Converse Scissors
- MMSU0238 - Adson Tissue Forceps, 12.5cm
- MMSU0241 - Adson Dissecting Forceps, 12.5cm

MMK043/2

Instruments





Malosa Medical

SINGLE-USE INSTRUMENTS

Cosmetic Surgery Packs

Liposuction



MMK045

Contents:

- MMSU3020 - BP No. 3 Scalpel Blade Handle with Gradations in cm
- MMSU0050 - Rampley Sponge Holding Forceps, 25cm
- MMSU0182 - Crile-Murray Needle Holder, 13cm
- MMSU0159 - Mayo Scissors, 16.5cm
- MMSU0238 - Adson Tissue Forceps, Lightweight, 12.5cm
- MMSU0218 - Standard Scissors, Blunt/Blunt, 13cm

Instruments



3020

BP No.3



0050

Rampley



0182

Crile-Murray



0159

Mayo



0238

Adson



0218

Standard



Ocular Surgery Packs

Minor Oculoplastics Surgery Packs

MMK085 Oculoplastics Pack

SELECT FROM OUR
STANDARD PACKS
OR
WE WILL CUSTOMISE
A PACK TO YOUR
OWN DESIGN



3001 - Desmarres Lid Retractor



1213 - Westcott Scissors



3012 - Castroviejo Needle Holder



0114/15 - #15 Blade Plastic Handled Scalpel



0276 - Bishops Forceps

3007 - Halstead Mosquito Forceps x 4

1254 - Jaffe Retractor



0068 - Plastic Holding Ramplery



0208 - Stephens Tenotomy Scissors

Chalazion Packs

Complete Chalazion Packs

MMK132 - CHALAZION PACK WITH DOUBLE ENDED CURETTE



1150 - Chalazion Clamp Small



1157 - Double Ended Curette 1 & 2mm




0114/11 - #11 Blade Disposable Scalpel



1149 - Chalazion Clamp - Medium

1148 - Chalazion Clamp - Large



1156 - 3.0mm Curette

1155 - 2.0mm Curette

1151 - 1.0mm Curette



Malosa Medical

SINGLE-USE INSTRUMENTS

Cosmetic Surgery Packs



Minor Plastics Pack

MMK0307

Contents:

- MMSU0249 Iris Dissecting Scissors, Straight, 11.5cm
- MMSU0183 Derf Needle Holder, without groove, 12cm
- MMSU0238 Adson Tissue Forceps, 1.25mm tip, 12.5cm
- MMSU0241 Adson Dissecting Forceps, 1.25mm tip, 12.5cm
- MMSU0114 #15 Blade Plastic Handled Scalpel
- MMSU3037 Gillies Skin Hook, Fine, 16.5cm x 2
- MMSU3007 Halstead Mosquito Forceps (Micro), Curved x 2
- MMSU0243 Curved Metzenbaum Scissors, 11.5cm
- MMSU0255 Mayo Scissors, 14.5cm
- MMSU0147 60ml Gallipot
- MMSU0120 Non Woven Swabs x 5
- MMSU0502 Round Fenestrated Drape
- MMSU0137 Supawrap 50 x 50
- MMSU0170 Plastic Procedure Tray




Minor Ops Hand Pack

MMK0310

Contents:

- MMSU0068 Plastic Holding Ramp, Straight, 18cm x 2
- MMSU0218 Standard Scissors, 13cm
- MMSU3037 Gillies Skin Hook, Fine, 16.5cm
- MMSU0182 Crille-Murray Needle Holder without groove, 13cm
- MMSU3020 Scalpel Blade Handle with Gradations in CM x 2
- MMSU0238 Adson Tissue Forceps, 1.25mm tip, 12.5cm
- MMSU0214 Stephens Tenotomy Scissors, Straight, Blunt tip
- MMSU0243 Curved Metzenbaum Scissors, 11.5cm
- MMSU3030 Lagenbach Retractor, 21.5cm handle
- MMSU3003 West Retractor, 3 x 4 Teeth, Sharp
- MMSU0170 Plastic Procedure Tray
- MMSU0137 Supawrap 50 x 50





Malosa Medical

SINGLE-USE INSTRUMENTS

Suture Packs

Suturing/Minor Ops



Contents:

- MMSU0141 - Curved Blade Scissors, Pointed 9cm
- MMSU0241 - Adson Dissecting Forceps, 12.5cm
- MMSU0238 - Adson Tissue Forceps, 12.5 cm
- MMSU0182 - Crile-Murray Needle Holder, 13cm
- MMSU0147 - 60ml Gallipot
- MMSU0120 - Non-woven Swabs, 10cm x 10cm
- MMSU0502 - Round Fenestrated Drape

MMK0304

Suture Removal



Contents:

- MMSU0141 - Curved Blade Scissors, Pointed 9cm
- MMSU0241 - Adson Dissecting Forceps, 12.5cm
- MMSU0147 - 60ml Gallipot
- MMSU0120 - Non-woven Swabs, 10cm x 10cm

MMK0303

Instruments

Instruments are placed on a plastic tray, wrapped in a 50 x 50 sterile field and pouched.

0141 9cm Long



Curved Blade Scissors

0182 13cm Long



Crile-Murray Needle Holder

0502



Fenestrated Drape

0238 & 0241 12.5cm Long



Adson Forceps

1.2mm Tip
1 x 2 Teeth



0238

1.2mm Tip



0241



Gallipot & Swabs



Malosa Medical

T: +44 (0)870 3000 555

F: +44 (0)870 3001 555

E : info@malosa.com

W: www.malosa.com