



# Abdominal Liposuction Cannula

Available in a range of tips, port configurations, lengths and diameters.



Pressure release vent on thumb grip

Bespoke designs available for a minimum order of 250pcs

### Pyramid



Ref #	Length (mm)	Dia (mm)	Port (mm)
3104	230	3.0	5.5 X 2.3
3106	260	4.0	7.0 X 2.5
3113	275	6.0	8.5 X 3.0

### Mercedes

#### Standard Tip



Ref #	Length (mm)	Dia (mm)	Port (mm)
3103	230	3.0	5.5 X 2.0
3105	260	4.0	7.0 X 2.25
3108	260	5.0	8.5 X 2.5

#### Bullet Tip



Ref #	Length (mm)	Dia (mm)	Port (mm)
3126	230	3.0	5.5 X 2.0
3125	260	4.0	7.0 X 2.0
3127	260	5.0	8.5 X 2.3



# Facial Liposuction Cannula

Cannula with leur lock hub, ports in a radial orientation, with a range of tips, port numbers, lengths and diameters.



### Standard Tip

#### Mono-Port




Ref #	Length (mm)	Dia (mm)	Port (mm)
3159	80	2.0	3.5 X 1.0

#### Bi-Port Radial



Ref #	Length (mm)	Dia (mm)	Port (mm)
3124	185	2.75	4.0 X 2.0

#### Tri-Port Radial



Ref #	Length (mm)	Dia (mm)	Port (mm)
3119	70	2.0	3.0 X 1.0
3120	110	3.0	4.0 X 1.5
3122	110	2.5	3.0 X 0.8
3123	185	2.75	4.0 X 2.0



# Facial Liposuction Cannula

Cannula with leur lock hub, ports aligned on the axis with a range of tips, port numbers, lengths and diameters.


### Standard Tip

#### Bi-Port Axial



Ref #	Length (mm)	Dia (mm)	Port (mm)
3179	150	2.0	3.0 X 1.5
3181	150	3.0	6.0 X 2.0

#### Tri-Port Axial



Ref #	Length (mm)	Dia (mm)	Port (mm)
3199	80	3.0	3.5 X 2.0


### Spatula Tip

#### Mono-Port



Ref #	Length (mm)	Dia (mm)	Port (mm)
3142	135	2.5	3.0 X 1.5

#### Bi-Port Axial



Ref #	Length (mm)	Dia (mm)	Port (mm)
3171	135	2.0	5.0 X 1.0
3172	185	2.0	5.0 X 1.0
3173	135	3.0	6.5 X 2.0
3174	185	3.0	6.5 X 2.0



### Liposuction Cannula

Cannula with leur lock for tumescent infiltration and fat repositioning

#### Fat Repositioning Cannula

Extra wide and longer port to prevent fat cell damage



Ref #	Length (mm)	Dia (mm)	Port (mm)
3225	80	2.0	4.0 X 1.5
3226	90	2.0	4.0 X 1.5

### Infiltration Cannula

#### 4 Spiral Port




Ref #	Length (mm)	Dia (mm)	Port (mm)
3218	200	3.0	1.0 Dia

#### 8 Spiral port



Ref #	Length (mm)	Dia (mm)	Port (mm)
3217	325	2.75	1.0 Dia
3219	185	2.0	1.0 Dia

### Leur Lock Connector



Ref #	Length (mm)	Dia (mm)	Port (mm)
3100	20	9.0	4.0 Dia



# Malosa Medical

## SINGLE-USE INSTRUMENTS

# Cosmetic Surgery Packs

## Liposuction



MMK045

**Contents:**

- MMSU3020 - BP No. 3 Scalpel Blade Handle with Gradations in cm
- MMSU0050 - Rampley Sponge Holding Forceps, 25cm
- MMSU0182 - Crile-Murray Needle Holder, 13cm
- MMSU0159 - Mayo Scissors, 16.5cm
- MMSU0238 - Adson Tissue Forceps, Lightweight, 12.5cm
- MMSU0218 - Standard Scissors, Blunt/Blunt, 13cm

## Instruments

