



Malosa Medical

SINGLE-USE INSTRUMENTS

Suture Packs

Suturing/Minor Ops



MMK0304

Contents:

- MMSU0141 - Curved Blade Scissors, Pointed 9cm
- MMSU0241 - Adson Dissecting Forceps, 12.5cm
- MMSU0238 - Adson Tissue Forceps, 12.5 cm
- MMSU0182 - Crile-Murray Needle Holder, 13cm
- MMSU0147 - 60ml Gallipot
- MMSU0120 - Non-woven Swabs, 10cm x 10cm
- MMSU0502 - Round Fenestrated Drape

Suture Removal



MMK0303

Contents:

- MMSU0141 - Curved Blade Scissors, Pointed 9cm
- MMSU0241 - Adson Dissecting Forceps, 12.5cm
- MMSU0147 - 60ml Gallipot
- MMSU0120 - Non-woven Swabs, 10cm x 10cm

Instruments

Instruments are placed on a plastic tray, wrapped in a 50 x 50 sterile field and pouched.

